

**BioPureSPN™ MACRO HIL-PSA Used as IEX Columns & 96-Well Plates**  
(50-200µL elution volume, 35-350 µg (HILIC) 0.2-1 mg (IEX) maximum capacity)

---

When used in an ion exchange mode, these primary-secondary amine spin columns, PSA™ ERLIC, will retain peptides and proteins with multiple negative charges and will remove non-ionic and basic detergents in the presence of 30% organic solvent. Polar and negatively charged solutes will be preferentially retained. Use for preliminary fractionation by negative and neutral charge differences for MS samples.



**Directions:** (p/n: **HMM HIL-PSA, HMM HIL-PSA.20, HNS HIL-PSA-L**): Cut Off The Outlet Tab with Side-Cut Pliers or a Box Cutter against an elevated support surface. Loosen The Cap when spinning.

- **Conditioning the column:** Pipette 400µL of conditioning solvent (e.g., 100% acetonitrile or MeOH) into the column and centrifuge for 1 min. at about 55 x g (@ ~400 rpm with an Eppendorf micro centrifuge). Flush with 1 tube volume (~400µL ea.) of 100% water. Flush with at least 2 bed volumes (200µL) of 100% water before elution with any salt solution to prevent salt precipitation. Condition the cartridge with a strong buffer for at least one hour prior to its initial use (add 200µL of conditioning buffer, spin for 10 sec. and let stand in the tube). A convenient solution to use is 0.2 M monosodium phosphate + 0.3 M sodium acetate (the pH will be between 3.0 and 6.5). Spin out remaining strong buffer and flush with 200µL, 100% water.
- **Equilibrating the column for IEX:** Pipette 200µL of 5mM phosphate (or 25mM ammonium formate), pH 6.5 and centrifuge it for 1 min. at about 30 x g (@~200 rpm with an Eppendorf micro centrifuge). Repeat two more times. Remove the collecting tube and blot dry any moisture on the exterior of the column.
- **Processing the sample:** (*Note: When using fixed-angle rotors, place a mark on the upper side of the column. Place column in centrifuge with the mark facing outward in all subsequent centrifugation steps. Improper orientation will result in reduced recovery efficiency.*) Load your 50-200µL of sample (in 5mM, pH 6.5 buffer) to the column and place it in a new 2ml centrifuge tube. Spin the tube 1 min. at 110x g. Peptides and proteins will be retained, while detergents and non-polar solutes will elute in the liquid in the collecting tube (if some organic solvent is present in the buffer). Discard (or save) this liquid. Rinse with 50-100µL of loading or equilibration buffer to wash out any traces of impurities from your sample of interest.
- **Releasing the sample:** Add 50-200µL of 5mM phosphate + 0.1 - 0.4M NaCl (or 500mM ammonium formate or some other volatile electrolyte) to completely wet the frit, preferably containing 5-20% ACN to reduce hydrophobic association. Spin as above. Peptides and proteins will be in the liquid in the collection tube. If a sample is especially ionic it may be necessary to repeat this step to elute all of the sample. If especially acidic (poly-phosphorylated peptides), repeat this step with increasing amounts of salt to elute all of the sample at pH 2.0.

**NOTES:**

- These spin columns of PSA™ ERLIC will retain anionic solutes such as phosphopeptides or proteins.

**BioPureSPN™ MACRO HIL-PSA Used as HILIC (ERLIC) Columns & 96-Well Plates**  
(50-200µL elution volume, 35-350 µg (HILIC) 0.2-1.0 mg (IEX) maximum capacity)

---

When used in a HILIC (ERLIC) mode, these primary-secondary amine spin columns, PSA™ ERLIC, will retain soluble peptides and will remove non-ionic and cationic detergents during the load in > 50% organic solvent. More non-polar, and positively charged solutes will be less retained than anionic ones. Use for preliminary fractionation by polarity differences for MS samples.



**Directions:** (p/n: *HMM HIL-PSA, HMM HIL-PSA.20, HNS HIL-PSA-L*): Cut Off The Outlet Tab with Side-Cut Pliers or a Box Cutter against an elevated support surface. Loosen The Cap when spinning.

- **Conditioning the column:** Pipette 400µL of conditioning solvent (e.g., 100% acetonitrile or MeOH) into the column and centrifuge for 1 min. at about 55 x g (@ ~400 rpm with an Eppendorf micro centrifuge).
- **Equilibrating the column for HILIC:** Flush with 1 tube volume of 85% ACN, 5mM - 20mM ammonium formate (or acetate) at a pH appropriate for the separation (e.g. pH 6.5 to enhance fractionation by neutral and negative charges, or pH 3 to retain components a neutral column surface.). Centrifuge for 1 min. at 30 x g (@ ~200 rpm on an Eppendorf micro centrifuge). Repeat twice. Remove the collecting tube and blot dry any moisture on the exterior of the column.
- **Processing the sample:** (*Note: When using fixed-angle rotors, place a mark on the upper side of the column. Place column in centrifuge with the mark facing outward in all subsequent centrifugation steps. Improper orientation will result in reduced recovery efficiency.*) Load 50-200µL of sample (in the same buffer at an appropriate pH) to the column and place it in a new 2mL centrifuge tube. Spin the tube 1 min. at 110x g. Polar analytes will be retained, while detergents and non-polar solutes will elute in the liquid in the collecting tube. Discard (or save) this liquid. Rinse with 50-100µL of loading or equilibration buffer to wash out any traces of impurities from your sample of interest.
- **Releasing the sample:** Add 50-200µL of 40% ACN, 20mM - 100mM ammonium formate or some other volatile electrolyte to completely wet the frit. Spin as above. Peptides and proteins will be in the liquid in the collection tube. If a sample is especially ionic it may be necessary to repeat this step to elute all of the sample. If especially acidic (poly-phosphorylated peptides), repeat this step with increasing amounts of salt to elute all of the sample at pH 2.0.

**NOTES:**

---

- For a discussion of the ERLIC technique see the ERLIC-WAX Dropbox on our web site: <http://www.nestgrp.com>

PSA™ ERLIC is a trademark of The Nest Group, Inc.  
12.03.18