

## Nucleotide Analysis with Vydac 302IC4.6 Ion-Chromatography Columns

Here's a unique bio research application for Vydac's 302IC4.6 column, originally designed for separating inorganic ions in environmental analysis. The packing, a low-capacity quaternary amine anion exchanger based on high-purity 10 $\mu$ m large-pore silica, separates organic ions by a mixture of ionic and reversed-phase (hydrophobic) interactions.

### Environmentally Friendly Mobile Phase: Easy Disposal

The nucleotide separation shown in Figure 1 was performed with a mobile phase that is atypical by reversed-phase standards because it contains no organic solvent, thus simplifying disposal. Instead, the acetic acid used for pH adjustment takes the place of the organic solvent in facilitating the reversed-phase portion of the mixed-mode separation.

### Retention and Mobile-Phase Effects

As indicated by the elution order, ionic interactions due to the negative charges of the phosphates appear to play a major role, separating mono-, di-, and tri-phosphate nucleotides as groups. Within each group, differing properties of the bases result in separation of individual nucleotides. Hydrophobic interactions vary with size of the heterocyclic ring structure and tend to cause purine nucleotides to be retained longer than pyrimidine nucleotides.

However, the bases also carry positive charges at the mobile-phase pH, thereby moderating ionic retention. Uracil is the least basic, reflecting its lack of an amino substituent on the heterocyclic ring. This accounts for longer retention of uridine nucleotides and their apparent anomalous elution between the corresponding purine nucleotides.

The 1:1 proportion of dibasic and monobasic phosphate salts in the mobile phase was important in achieving conditions for this separation. Using only dibasic phosphate eluted all peaks too early, probably because more acetic acid was needed to adjust the pH. Using exclusively monobasic phosphate eluted all peaks much later.








Figure 1. Separation of twelve common nucleotides on Vydac ion chromatography column. Column: Vydac 302IC4.6, 4.6mmID x 250mmL. Flow rate: 2.0 mL/min. Detection: 260nm. Mobile phase: A = NaH<sub>2</sub>PO<sub>4</sub>/Na<sub>2</sub>HPO<sub>4</sub> (1:1 molar ratio) at 25mM total concentration in water, adjusted to pH 2.8 with acetic acid. B = NaH<sub>2</sub>PO<sub>4</sub>/Na<sub>2</sub>HPO<sub>4</sub> (1:1 molar ratio) at 125mM total concentration in water, adjusted to pH 2.9 with acetic acid. Gradient: 0% B for 2 minutes, then linear from 0% to 100% B in 17 minutes, hold 100% for 2 minutes, return to 0% B in 0.1 minute. Sample: 3 $\mu$ L containing 0.25 $\mu$ g of each nucleotide.

#### pyrimidine nucleotides



#### purine nucleotides



### Ordering Information

Cat. No.	Description
302IC4.6	Column, Ion Chromatography, Quarternary Amine, 10 $\mu$ m, 4.6mm ID x 250mm L

To place an order, call The Nest Group, 800.347.6378 your local Vydac distributor.

Copyright © 1998, The Separations Group, Hesperia, CA, USA. All rights reserved. "Vydac" is a trademark of The Separations Group. Prices and specifications subject to change without notice.

*an ISO 9001 Certified company*

*serving chromatographers with quality  
packings and bonding chemistry since 1971*

The Nest Group • 45 Valley Road • Southborough, MA 01772.1306 • USA  
Toll-Free: 800.347.6378 • Tel: 508.481.6223 • Fax: 508.485.5736 • EMail: info@nestgrp.com