


A photograph of a field filled with numerous haystacks. The haystacks are made of dry, golden-brown straw and are secured with wooden stakes. The field is green with grass, and the background shows rolling hills and mountains under a blue sky with some clouds. The image is framed by a white border with small black corner marks.

COMBINATORIAL LIBRARY SCREENING

**Discover Big
Opportunities
In Small
Molecules**

Drug-Protein Conjugate Fast SEC Columns:



CONDITIONS

Column: PolyHYDROXYETHYL A,
2.1 x 35 mm

Mobile Phase: 50 mM NaPO₄, pH 7.2
with 2% DMSO (v/v)

Flow Rate: 0.4 mL/min

Sample: 2 μ L with 50 μ M
of each solute

Data courtesy of Krishna Kalghargi,
NewGenesis Inc. (Cambridge, MA)

- Isolate conjugates in about 20 seconds.
- Use as little as 10 nM of sample.
- Separate conjugates from small solutes with a size exclusion limit of <1 KDa.
- Non-interactive even with weak binding constants.
- Couple to reverse phase columns for LC-MS.

ORDERING INFORMATION

Part Number	Description	Price
P051HY05006	1.0 x 50 mm, 5 μ m	\$445.00
P03.52HY05006	2.1 x 35 mm, 5 μ m	\$445.00
P03.54HY05006	4.6 x 35 mm, 5 μ m	\$445.00

To place an order or for more information, call **1-800-347-6378**, e-mail info@nestgrp.com visit <http://world.std.com/~nestgrp>, or add your phone number and e-mail address to the space provided and **fax this side of the card to 1-508- 485-5736**.

The Nest Group

The Nest Group Inc.
45 Valley Road
Southborough, MA 01772

TELEPHONE _____

E-MAIL _____