

Glycoprofiling Strategies to Support Production of Recombinant Therapeutic Proteins

Shekar Ganesa, Ph.D.
Associate Director
BioAnalytical Development
Genzyme Corporation
Framingham, MA 01701

The Genzyme logo is displayed in a green, lowercase, sans-serif font. The letters are closely spaced, and the 'g' and 'z' have a distinctive shape. The logo is positioned in the bottom right corner of the slide.

genzyme

Why Monitor Glycosylation?

- Impacts bioavailability/clearance due to terminal sialic acid and galactose residues
- Effects folding (FSH, LH, CG), protease susceptibility (Fibronectin), aggregation (Fibrinogen) and stability (Thyroxine-binding globulin) (Varki, 1993)
- Influences biological activity (EPO, FSH etc.)
- Ensures product batch-to-batch consistency and thus clinical efficacy

Glycosylation in Antibodies

Lund et. al., 1996 and Jefferies et. al., 1998


- Influences conformation and thus affects folding
- GlcNAc₂-Man important in Fc receptor binding and is involved in complement-mediated lysis


Davies et. al., 2001 and Shields et. al., 2002

- Lower fucose and bisecting GlcNAc content improves FcγRIIIA receptor binding and enhances antibody-dependent cell-mediated cytotoxicity (ADCC)


Kaneko et. al., 2006

- G0 variant is pro-inflammatory and Fc sialylation is anti-inflammatory





Fucose may hinder the approach of glycans on Asn-162 of FcγRIII, thereby reducing interaction between IgG and the FcγRIII (Ferrara et. al., 2006)




Depletion of Fucose from IgG Fc may allow Tyr-296 to interact with Lys-128 in FcγRIII, thereby increasing affinity (Okazaki et. al., 2004)


Factors Influencing Glycosylation

- Expression system (CHO, NS0, Pichia, Plants, Baculovirus)
- Cell age
- pH
- Cell density
- pCO₂
- Ammonium concentration and other metabolites (i.e. sodium butyrate, Jenkins et. al., 1996)
- Downstream processing


Expression System and Glycosylation



CHO




Baculovirus




CHO + Kifunensine

vanPatten et al. *Glycobiology*, 2007 (accepted)



CHO - Lec1 mutant *genzyme*

Molecular Complexity - Example



~ 4.0 x 10⁴ possible glycovariants

- Robust analytical platforms are critical to assess product quality and consistency

Monitoring Glycosylation

- Individual glycan structure
- Intact protein glycoforms
- Site-specific glycosylation


Monitoring Individual Glycan Structure

- Monosaccharide compositional analyses
 - Monosaccharides (Fuc, Gal, GlcNAc, Man)
 - Sialic acid
 - Mannose-6-phosphate
- Oligosaccharide Profiling
 - Oligosaccharide profiling monitors potential changes in oligosaccharide structures and their relative population

Oligosaccharide Profiling Methods


- Unlabeled glycans
 - High pH Anion exchange chromatography with pulsed amperometric detection (HPAEC-PAD)
 - Mass Spectrometry (MALDI-TOF)
- Fluorescently labeled glycans
 - Anthranilic acid, AA \longrightarrow NP-HPLC
 - 2-aminobenzoic acid, 2-AB \longrightarrow NP-HPLC
 - 8-aminopyrene-1,3,6-trisulfonate, APTS \longrightarrow CE-LIF

Unlabeled Glycans by HPAEC-PAD




Column : Dionex PA-100 column (250 X 4.6 mm) at 35 °C
Mobile Phase A: 100 mM Sodium Hydroxide
Mobile Phase B: 100 mM Sodium Hydroxide + 500 mM Sodium Acetate
Gradient: 0%-55%B in 85 minutes
Detection: Pulsed Amperometric Detection (4-stage waveform)



Unlabeled Glycans by MALDI-TOF




mAb: MALDI-TOF of Unlabeled Glycans




Fluorescence Labeling of Glycans



AA



2-AB



APTS


genzyme

*Anumula et al. *Glycobiology* 1998, 8, 685

*S. Ma, *Anal. Chem.* 1999, 71, 5185


*Bigge et al. *Anal. Biochem.* 1995, 230, 229

Oligosaccharide Profiling by AA-labeling/NP-HPLC




Column : Luna NH2, 5 μm , 100 \AA (250 x 4.6 mm) at 50 $^{\circ}\text{C}$
Mobile Phase A: 2% acetic acid, 1% THF in acetonitrile
Mobile Phase B: 5% acetic acid, 3% TEA, 1% THF in water
Gradient: 30%-95%B in 85 minutes
Detection: ex. 230 nm and em. at 425 nm

mAb: AA-labeling/NP-HPLC




Column : Asahipak NH2P-50 4D (150 x 4.6 mm) at 50 °C
Mobile Phase A: 2% acetic acid, 1% THF in acetonitrile
Mobile Phase B: 5% acetic acid, 3% TEA, 1% THF in water
Gradient: 30%-50%B in 60 minutes
Detection: ex. 230 nm and em. at 425 nm

mAb: 2-AB labeling/NP-HPLC



Column : Polysulfoethyl A (100 x 4.6 mm) at 45 °C
Mobile Phase A: Acetonitrile
Mobile Phase B: Water
Gradient: 15%-41%B in 42 minutes
Detection: ex. 330 nm, em. 420 nm

mAb: APTS-labeling/CE-LIF



Column : eCAP™ N-CHO capillary (50 μm x 27 cm) at 20 °C
Cathode buffer: 2% acetic acid, 1% THF in acetonitrile
Anode buffer: 5% acetic acid, 3% TEA, 1% THF in water
Conditions: Reverse polarity, 20 kV, 20 μA
Detection: ex. 488 nm and em. at 520 nm


HPAEC-PAD versus Fluorescence

	<u>HPAEC-PAD</u>	<u>Fluorescence</u>
Charge/linkage based rugged separation	Yes	Yes
Sample derivatization and clean-up (prior to analysis)	No	Yes
Quantitative	No	Yes
Assay variability	High	Low
MS compatible	Maybe	Yes
		enzyme


Monitoring Intact Protein Glycosylation

- Isoelectric focusing
 - Snapshot to examine charge heterogeneity
 - No structural detail obtained
- Mass Spectrometry (i.e. for mAbs)
 - Powerful characterization tool that provides a snapshot of glycosylation at molecular level
 - Eliminates potential selectivity associated with glycan release
 - Quantitative based on the assumption of equal ionization efficiencies for various glycoforms

Intact Protein Glycosylation : IEF Analysis




- pI shifts with removal of sialic acid




- Sialylation affects Specific Activity of this protein


Intact Protein Glycosylation (mAb) – Mass Spectrometry




ESI-MS Spectra of mAb Heavy Chain



Monitoring Site-Specific Glycosylation




Trpytic Map of a Multiply Glycosylated Protein



High Degree of Lot to Lot Consistency





- Relative intensities of various glycans at Asn-X in 7 Drug Substance lots







Site-Specific Glycosylation (LC/MS)

- Advantages
 - Detailed site-specific structural analysis and relative abundance of different glycan species at each site
 - Semi-quantitative and eliminates potential selectivity associated with glycan release
- Disadvantages
 - Assumes ionization efficiencies for all glycopeptides similar
 - Requires high level of technical expertise and instruments (ion trap, Q-TOF etc.)

Precision : HPAEC-PAD versus AA-labeling


HPAEC-PAD	% Peak Area (N=39)				
		+			
Average	52.3			21.3	19.6
SD	3.3			3.2	3.0
%CV	12.6			15.0	15.3

AA/HPLC	% Peak Area (N=30)				
		+			
Average	50.2			24.4	20.9
SD	1.0			1.1	1.3
%CV	2.1			4.5	6.4


genzyme

APTS-labeling/CE-LIF versus MALDI-TOF

- Relative percentage of G0 glycans



AA versus MALDI-TOF versus Intact MW



$$\text{Percent Galactosylation} = \frac{\text{Total Galactose Residues}}{\text{Total sites for Galactosylation}}$$

Conclusions

- Oligosaccharide profiling by HPAEC-PAD has significant limitations i.e. higher variability, non-quantitative, lack of structural information etc.
- Fluorescence-based methods offer a rugged, quantitative, alternative to HPAEC-PAD with on-line detailed structural information (LC-MS and MS/MS)
- Monosaccharides critical for initial characterization of the oligosaccharide structures must be routinely monitored
- The ultimate choice of any method depends on the glycoprotein being analyzed, analyst expertise and the successful implementation of that particular technology

Acknowledgments

John Harrahy

Kelly Toler

Brad Whittaker

Rebecca Sendak

Kefei Zheng

Jennifer Baker-Malcolm

Anders Lund

Daotian Fu

Kate Zhang

Xiaoying Jin

Qun Zhou

Scott vanPatten

Karen Lee

Tim Edmunds

Bob Mattaliano

Joseph Siemiatkoski, Biogen-Idec

Stacey Ma, Genentech