


# Combinatorial Library

## *FAST* Protein Conjugate Screening

Reduces Protein Requirements 100 X

**20 Second Assays!**



*Data courtesy of Krishna Kalghatgi, NeoGenesis Inc. (Cambridge, MA)*

PolyHYDROXYETHYL A™, 4.6 x 50mm, 5µm, 60Å  
50 mM Na-PO<sub>4</sub>, Ph 7.2, with 2% DMSO (v/v) Mobile Phase  
2 mL/min Flow Rate  
A<sub>254</sub> Detection  
10µL Sample with 50 µM of each solute