

# Fast Macro DispoDIALYZER™

Samples from 1 ml to 10 ml (Single use)

## advantages

- Macro-volume dialysis
- Economical
- Easy to use
- Leak-proof
- Low protein binding
- High sample recovery
- Rapid dialysis/purification

## applications

- Exchange of buffers
- Removal of detergents
- Concentration of samples
- HPLC, HPCE
- Purification of proteins, DNA and RNA
- Removal of excess radiolabel or PCR-primers
- GC, GC-MS, NMR



Fast Macro DispoDIALYZER™ Assembled

### Fast Macro DispoDIALYZER



The Fast Macro DispoDIALYZER is also available for larger samples. This unit is available with 5,000 or 10,000 Dalton regenerated cellulose membranes. The other basic features of this unit are the same as those of the smaller volume DispoDIALYZERS described earlier. The Fast Macro DispoDIALYZER includes shell, two rings and two membranes of MWCO 5KDa or 10KDa (4.5 cm<sup>2</sup> surface area each). Assembly required.

## Fast Macro DispoDIALYZERS – Ordering Information

Membrane MWCO (Daltons)	Qty. of 25	Qty. of 50	Qty. of 100
5,000	SDIS 0525K	SDIS 0550K	SDIS 0510K
10,000	SDIS 1025K	SDIS 1050K	SDIS 1010K