

CAPELLINI 6-32

CAPILLARY HPLC COLUMNS with 6-32 Thread PEEK/Ti Fittings



Capillary HPLC Columns

Unique Features:

- 6-32 threads for capillary fittings
- PEEK/Titanium fittings
- Tested and supplied with CofA
- High performance
- Higgins Analytical Sorbents

Applications:

- LC-MS
- Microbore HPLC

Standard Sizes:

- 50 to 150 x 0.150mm
- 50 to 150 x 0.075mm
- Longer and shorter lengths available

Guide to CAPELLINI Part Numbers

xFS-xx07-xxxx 75µm CAPELLINI 6-32
 xFS-xx15-xxxx 150µm CAPELLINI 6-32

See Page 23 for complete Part Number information

Higgins Analytical Capellini are robust, high performance capillary HPLC columns for LC-MS and other very small scale applications. The inert PEEK and titanium fittings complement the Higgins Analytical sorbents such as TARGA which was developed to enable LC-MS applications without TFA.


Our no nonsense testing procedure does not resort to sample and flow splitting techniques, inaccurate k' values, and miniscule chromatograms, as do some manufacturers, to illustrate column performance. The CofA's symmetrical peaks and the high resolution map speak for themselves.

Capellini - Back by Popular Demand

When Higgins Analytical introduced Capillary columns three years ago, we removed Capellini from our catalog, however, we continued to manufacture this design for customers that had been using it for several years. Somewhat to our surprise, Capellini sales have continued to grow and we still receive numerous product inquiries in spite of their absence from our literature. While Capillary columns are very efficient and cost effective (see page 40), we have found that many chromatographers prefer the symmetrical, flow reversible design and the 6-32 threaded end-fittings of Capellini. Even though Capellini are more expensive to manufacture and somewhat less efficient than Capillary columns, due to the end fitting design, they enjoy a wide variety of uses that involve flow reversal, such trap applications.

High Selectivity:

A partial 1.5µL/min LCMSMS map of aBSA tryptic digest on a 150cm x 75µm Targa C18 5µm Capellini. P/N TC-1507-C185



CLIQUEUS Capillary Columns with 6-32 threads

| | | | | | | |
|------------------------|---|-------|---------------|-------|---------------|-------|
| CLIQUEUS C8 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | CFS-0507-C085 | \$435 | CFS-1007-C085 | \$435 | CFS-1507-C085 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | CFS-0515-C085 | \$435 | CFS-1015-C085 | \$435 | CFS-1515-C085 | \$435 |

| | | | | | | |
|-------------------------|---|-------|---------------|-------|---------------|-------|
| CLIQUEUS C18 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | CFS-0507-C185 | \$435 | CFS-1007-C185 | \$435 | CFS-1507-C185 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | CFS-0515-C185 | \$435 | CFS-1015-C185 | \$435 | CFS-1515-C185 | \$435 |

| | | | | | | |
|---------------------------|---|-------|---------------|-------|---------------|-------|
| CLIQUEUS Cyano 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | CFS-0507-CNP5 | \$435 | CFS-1007-CNP5 | \$435 | CFS-1507-CNP5 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | CFS-0515-CNP5 | \$435 | CFS-1015-CNP5 | \$435 | CFS-1515-CNP5 | \$435 |

| | | | | | | |
|----------------------------|---|-------|---------------|-------|---------------|-------|
| CLIQUEUS Phenyl 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | CFS-0507-PHN5 | \$435 | CFS-1007-PHN5 | \$435 | CFS-1507-PHN5 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | CFS-0515-PHN5 | \$435 | CFS-1015-PHN5 | \$435 | CFS-1515-PHN5 | \$435 |

| | | | | | | |
|----------------------------|---|-------|---------------|-------|---------------|-------|
| CLIQUEUS Silica 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | CFS-0507-SIL5 | \$435 | CFS-1007-SIL5 | \$435 | CFS-1507-SIL5 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | CFS-0515-SIL5 | \$435 | CFS-1015-SIL5 | \$435 | CFS-1515-SIL5 | \$435 |

HAISIL100 Capillary Columns with 6-32 threads

| | | | | | | |
|--------------------------|---|-------|---------------|-------|---------------|-------|
| HAISIL 100 C8 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | HFS-0507-M085 | \$435 | HFS-1007-M085 | \$435 | HFS-1507-M085 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | HFS-0515-M085 | \$435 | HFS-1015-M085 | \$435 | HFS-1515-M085 | \$435 |

| | | | | | | |
|---------------------------|---|-------|---------------|-------|---------------|-------|
| HAISIL 100 C18 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | HFS-0507-M185 | \$435 | HFS-1007-M185 | \$435 | HFS-1507-M185 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | HFS-0515-M185 | \$435 | HFS-1015-M185 | \$435 | HFS-1515-M185 | \$435 |

| | | | | | | |
|------------------------------|---|-------|---------------|-------|---------------|-------|
| HAISIL 100 Silica 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | HFS-0507-SIL5 | \$435 | HFS-1007-SIL5 | \$435 | HFS-1507-SIL5 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | HFS-0515-SIL5 | \$435 | HFS-1015-SIL5 | \$435 | HFS-1515-SIL5 | \$435 |



CAPELLINI 6-32 (cont.)

CAPILLARY HPLC COLUMNS with 6-32 Thread PEEK/Ti Fittings



HAISIL300 Capillary Columns with 6-32 threads

| HAISIL 300 C18 5µm | 75µm and 150µm Capillary columns | | | | | |
|--------------------|----------------------------------|-------|---------------|-------|---------------|-------|
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | HFS-0507-W185 | \$435 | HFS-1007-W185 | \$435 | HFS-1507-W185 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | HFS-0515-W185 | \$435 | HFS-1015-W185 | \$435 | HFS-1515-W185 | \$435 |

HAISIL HL Capillary Columns with 6-32 threads

| HAISIL HL C18 5µm | 75µm and 150µm Capillary columns | | | | | |
|-------------------|----------------------------------|-------|---------------|-------|---------------|-------|
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | LFS-0507-C185 | \$435 | LFS-1007-C185 | \$435 | LFS-1507-C185 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | LFS-0515-C185 | \$435 | LFS-1015-C185 | \$435 | LFS-1515-C185 | \$435 |

HEAVY Capillary Columns with 6-32 threads

| HEAVY C18 5µm | 75µm and 150µm Capillary columns | | | | | |
|---------------|----------------------------------|-------|---------------|-------|---------------|-------|
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | VFS-0507-C185 | \$435 | VFS-1007-C185 | \$435 | VFS-1507-C185 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | VFS-0515-C185 | \$435 | VFS-1015-C185 | \$435 | VFS-1515-C185 | \$435 |

PROTO 200 Capillary Columns with 6-32 threads

| PROTO 200 C4 5µm | 75µm and 150µm Capillary columns | | | | | |
|------------------|----------------------------------|-------|---------------|-------|---------------|-------|
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | RFS-0507-D045 | \$435 | RFS-1007-D045 | \$435 | RFS-1507-D045 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | RFS-0515-D045 | \$435 | RFS-1015-D045 | \$435 | RFS-1515-D045 | \$435 |

| PROTO 200 C18 5µm | 75µm and 150µm Capillary columns | | | | | |
|-------------------|----------------------------------|-------|---------------|-------|---------------|-------|
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | RFS-0507-D185 | \$435 | RFS-1007-D185 | \$435 | RFS-1507-D185 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | RFS-0515-D185 | \$435 | RFS-1015-D185 | \$435 | RFS-1515-D185 | \$435 |



PROTO 300 Capillary Columns with 6-32 threads

| | | | | | | |
|-------------------------|---|---------------|---------------|---------------|---------------|-------|
| PROTO 300 C4 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | RFS-0507-W045 | \$435 | RFS-1007-W045 | \$435 | RFS-1507-W045 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| RFS-0515-W045 | \$435 | RFS-1015-W045 | \$435 | RFS-1515-W045 | \$435 | |

| | | | | | | |
|--------------------------|---|---------------|---------------|---------------|---------------|-------|
| PROTO 300 C18 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | RFS-0507-W185 | \$435 | RFS-1007-W185 | \$435 | RFS-1507-W185 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| RFS-0515-W185 | \$435 | RFS-1015-W185 | \$435 | RFS-1515-W185 | \$435 | |

PHALANX Capillary Columns with 6-32 threads

| | | | | | | |
|------------------------|---|---------------|---------------|---------------|---------------|-------|
| PHALANX C18 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | PFS-0507-C185 | \$435 | PFS-1007-C185 | \$435 | PFS-1507-C185 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| PFS-0515-C185 | \$435 | PFS-1015-C185 | \$435 | PFS-1515-C185 | \$435 | |

TARGA Capillary Columns with 6-32 threads

| | | | | | | |
|---------------------|---|---------------|---------------|---------------|---------------|-------|
| TARGA C8 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | TFS-0507-C085 | \$435 | TFS-1007-C085 | \$435 | TFS-1507-C085 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| TFS-0515-C085 | \$435 | TFS-1015-C085 | \$435 | TFS-1515-C085 | \$435 | |

| | | | | | | |
|----------------------|---|---------------|---------------|---------------|---------------|-------|
| TARGA C18 3µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | TFS-0507-C183 | \$435 | TFS-1007-C183 | \$435 | TFS-1507-C183 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| TFS-0515-C183 | \$435 | TFS-1015-C183 | \$435 | TFS-1515-C183 | \$435 | |

| | | | | | | |
|----------------------|---|---------------|---------------|---------------|---------------|-------|
| TARGA C18 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | TFS-0507-C185 | \$435 | TFS-1007-C185 | \$435 | TFS-1507-C185 | \$435 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| TFS-0515-C185 | \$435 | TFS-1015-C185 | \$435 | TFS-1515-C185 | \$435 | |



CAPELLINI 10-32

CAPILLARY HPLC COLUMNS with 10-32 Thread s.s. Fittings



Capillary HPLC Columns

Unique Features:

- 10-32 threads for standard fittings
- s.s. End Fittings/Titanium Frits
- Tested and supplied with CofA
- High performance
- Higgins Analytical Sorbents

Applications:

- LC-MS
- Microbore HPLC

Standard Sizes:

- 50 to 150 x 0.150mm
- 50 to 150 x 0.075mm
- Longer and shorter lengths available

NEW

Higgins Analytical Capellini have been popular with chromatographers requiring an easy-to-use and robust capillary HPLC column. The original Capellini (now called Capellini 6-32) were available only with 6-32 fittings. Now, we have another version for chromatographers preferring 10-32 thread fittings - Capellini 10-32. Another difference between the two versions of Capellini is that the 10-32 threads are stainless steel rather than PEEK. Both types are made with titanium frits on each end, providing inertness and flow direction reversability.


Guide to CAPELLINI 10-32 Part Numbers

xPS-xx07-xxxx 75µm CAPELLINI 10-32
 xPS-xx15-xxxx 150µm CAPELLINI 10-32

See Page 23 for complete Part Number information

Robust design with High Performance:

To the right, aBSA tryptic digest on a 150cm x 75µm Targa C18 5µm Capellini. P/N TC-1507-C185



CLIQUEUS Capillary Columns with 10-32 threads

| | | | | | | |
|------------------------|---|---------------|---------------|---------------|---------------|-------|
| CLIQUEUS C8 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | CPS-0507-C085 | \$450 | CPS-1007-C085 | \$450 | CPS-1507-C085 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| CPS-0515-C085 | \$450 | CPS-1015-C085 | \$450 | CPS-1515-C085 | \$450 | |

| | | | | | | |
|-------------------------|---|---------------|---------------|---------------|---------------|-------|
| CLIQUEUS C18 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | CPS-0507-C185 | \$450 | CPS-1007-C185 | \$450 | CPS-1507-C185 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| CPS-0515-C185 | \$450 | CPS-1015-C185 | \$450 | CPS-1515-C185 | \$450 | |

| | | | | | | |
|---------------------------|---|---------------|---------------|---------------|---------------|-------|
| CLIQUEUS Cyano 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | CPS-0507-CNP5 | \$450 | CPS-1007-CNP5 | \$450 | CPS-1507-CNP5 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| CPS-0515-CNP5 | \$450 | CPS-1015-CNP5 | \$450 | CPS-1515-CNP5 | \$450 | |

| | | | | | | |
|----------------------------|---|---------------|---------------|---------------|---------------|-------|
| CLIQUEUS Phenyl 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | CPS-0507-PHN5 | \$450 | CPS-1007-PHN5 | \$450 | CPS-1507-PHN5 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| CPS-0515-PHN5 | \$450 | CPS-1015-PHN5 | \$450 | CPS-1515-PHN5 | \$450 | |

| | | | | | | |
|----------------------------|---|---------------|---------------|---------------|---------------|-------|
| CLIQUEUS Silica 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | CPS-0507-SIL5 | \$450 | CPS-1007-SIL5 | \$450 | CPS-1507-SIL5 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| CPS-0515-SIL5 | \$450 | CPS-1015-SIL5 | \$450 | CPS-1515-SIL5 | \$450 | |

HAISIL100 Capillary Columns with 10-32 threads

| | | | | | | |
|--------------------------|---|---------------|---------------|---------------|---------------|-------|
| HAISIL 100 C8 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | HPS-0507-M085 | \$450 | HPS-1007-M085 | \$450 | HPS-1507-M085 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| HPS-0515-M085 | \$450 | HPS-1015-M085 | \$450 | HPS-1515-M085 | \$450 | |

| | | | | | | |
|---------------------------|---|---------------|---------------|---------------|---------------|-------|
| HAISIL 100 C18 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | HPS-0507-M185 | \$450 | HPS-1007-M185 | \$450 | HPS-1507-M185 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| HPS-0515-M185 | \$450 | HPS-1015-M185 | \$450 | HPS-1515-M185 | \$450 | |

| | | | | | | |
|------------------------------|---|---------------|---------------|---------------|---------------|-------|
| HAISIL 100 Silica 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | HPS-0507-SIL5 | \$450 | HPS-1007-SIL5 | \$450 | HPS-1507-SIL5 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| HPS-0515-SIL5 | \$450 | HPS-1015-SIL5 | \$450 | HPS-1515-SIL5 | \$450 | |



CAPELLINI 10-32 (cont.)

CAPILLARY HPLC COLUMNS with 10-32 Thread s.s. Fittings



HAISIL300 Capillary Columns with 10-32 threads

| HAISIL 300 C18 5µm | 75µm and 150µm Capillary columns | | | | | |
|--------------------|----------------------------------|-------|---------------|-------|---------------|-------|
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | HPS-0507-W185 | \$450 | HPS-1007-W185 | \$450 | HPS-1507-W185 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | HPS-0515-W185 | \$450 | HPS-1015-W185 | \$450 | HPS-1515-W185 | \$450 |

HAISIL HL Capillary Columns with 10-32 threads

| HAISIL HL C18 5µm | 75µm and 150µm Capillary columns | | | | | |
|-------------------|----------------------------------|-------|---------------|-------|---------------|-------|
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | LPS-0507-C185 | \$450 | LPS-1007-C185 | \$450 | LPS-1507-C185 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | LPS-0515-C185 | \$450 | LPS-1015-C185 | \$450 | LPS-1515-C185 | \$450 |

HEAVY Capillary Columns with 10-32 threads

| HEAVY C18 5µm | 75µm and 150µm Capillary columns | | | | | |
|---------------|----------------------------------|-------|---------------|-------|---------------|-------|
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | VPS-0507-C185 | \$450 | VPS-1007-C185 | \$450 | VPS-1507-C185 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | VPS-0515-C185 | \$450 | VPS-1015-C185 | \$450 | VPS-1515-C185 | \$450 |

PROTO 200 Capillary Columns with 10-32 threads

| PROTO 200 C4 5µm | 75µm and 150µm Capillary columns | | | | | |
|------------------|----------------------------------|-------|---------------|-------|---------------|-------|
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | RPS-0507-D045 | \$450 | RPS-1007-D045 | \$450 | RPS-1507-D045 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | RPS-0515-D045 | \$450 | RPS-1015-D045 | \$450 | RPS-1515-D045 | \$450 |

| PROTO 200 C18 5µm | 75µm and 150µm Capillary columns | | | | | |
|-------------------|----------------------------------|-------|---------------|-------|---------------|-------|
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | RPS-0507-D185 | \$450 | RPS-1007-D185 | \$450 | RPS-1507-D185 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | RPS-0515-D185 | \$450 | RPS-1015-D185 | \$450 | RPS-1515-D185 | \$450 |



NEW

PROTO 300 Capillary Columns with 10-32 threads

| | | | | | | |
|-------------------------|---|-------|---------------|-------|---------------|-------|
| PROTO 300 C4 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | RPS-0507-W045 | \$450 | RPS-1007-W045 | \$450 | RPS-1507-W045 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | RPS-0515-W045 | \$450 | RPS-1015-W045 | \$450 | RPS-1515-W045 | \$450 |

| | | | | | | |
|--------------------------|---|-------|---------------|-------|---------------|-------|
| PROTO 300 C18 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | RPS-0507-W185 | \$450 | RPS-1007-W185 | \$450 | RPS-1507-W185 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | RPS-0515-W185 | \$450 | RPS-1015-W185 | \$450 | RPS-1515-W185 | \$450 |

PHALANX Capillary Columns with 10-32 threads

| | | | | | | |
|------------------------|---|-------|---------------|-------|---------------|-------|
| PHALANX C18 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | PPS-0507-C185 | \$450 | PPS-1007-C185 | \$450 | PPS-1507-C185 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | PPS-0515-C185 | \$450 | PPS-1015-C185 | \$450 | PPS-1515-C185 | \$450 |

TARGA Capillary Columns with 10-32 threads

| | | | | | | |
|---------------------|---|-------|---------------|-------|---------------|-------|
| TARGA C8 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | TPS-0507-C085 | \$450 | TPS-1007-C085 | \$450 | TPS-1507-C085 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | TPS-0515-C085 | \$450 | TPS-1015-C085 | \$450 | TPS-1515-C085 | \$450 |

| | | | | | | |
|----------------------|---|-------|---------------|-------|---------------|-------|
| TARGA C18 3µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | TPS-0507-C183 | \$450 | TPS-1007-C183 | \$450 | TPS-1507-C183 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | TPS-0515-C183 | \$450 | TPS-1015-C183 | \$450 | TPS-1515-C183 | \$450 |

| | | | | | | |
|----------------------|---|-------|---------------|-------|---------------|-------|
| TARGA C18 5µm | 75µm and 150µm Capillary columns | | | | | |
| | 50x0.075mm | | 100x0.075mm | | 150x0.075mm | |
| | TPS-0507-C185 | \$450 | TPS-1007-C185 | \$450 | TPS-1507-C185 | \$450 |
| | 50x0.150mm | | 100x0.150mm | | 150x0.150mm | |
| | TPS-0515-C185 | \$450 | TPS-1015-C185 | \$450 | TPS-1515-C185 | \$450 |

