


Can Your Metabolomics Column Retain Polar Analytes?


TARGA C18 Increases Resolution of Hydrophilic Analytes

Better Retention




Column: TARGA C18, p/n: H TS-2546-C185
Solvent: 100% water. Detection: LC-MS

Increased Resolution



Column: TARGA C18, p/n: H TS-2546-C185
Gradient: 100% water, 0.1% FA to 60% ACN

Better Capacity



Trapping Data Courtesy of Brian Hampton, Univ. Maryland

TARGA Piccolo μ Trap™ C18, 35 x 0.3mm
Loading at 98% water, 0.1% FA 2% ACN

Load Samples in 100% water
Retain Polar Analytes for Increased Resolution
10x Faster Equilibration in 100% Water for Shorter Cycle Times