


Combinatorial Library
***FAST* Protein Conjugate Screening**

Reduces Protein Requirements 100 X

20 Second Assays!



Data courtesy of Krishna Kalghatgi, NeoGenesis Inc. (Cambridge, MA)

PolyHYDROXYETHYL A™, 4.6 x 50mm, 5µm, 60Å
50 mM Na-PO₄, pH 7.2, with 2% DMSO (v/v) Mobile Phase
2 mL/min Flow Rate
A₂₅₄ Detection
10µL Sample with 50 µM of each solute