





iHILIC[®]

Advancing HILIC Separations at High pH



- Zwitterionic charge modulated amide/diol polymeric HILIC columns
- Complementary selectivities for separation of polar compounds
- Excellent durability and ultra-low column bleeding
- Universal for LC-MS based "Omics" studies and other applications at pH 1-10
- 5 μm in PEEK, SS, and PEEK-SS hardware



HILICON AB

Email: info@hilicon.com | Website: www.hilicon.com

©2023 HILICON AB. All rights reserved. | iHILIC[®] is a registered trademark of HILICON AB, Sweden

