

# The Nest Group, Inc.

## High Speed Separation of Catecholamines

### Maccel 120-3-C18 AQ

Part number: B0404F184PS030  
Dimension: 33 × 4 mm  
Eluent: CAT-A-Phase  
Flow: 1.0 ml/min  
Detection: Coulometric  
Temperature: 30 °C  
Injection: 5 µl  
Sample:  
1: Norepinephrine  
2: Epinephrine  
3: Dihydroxybenzylamine  
4: Dopamine

### Maccel 120-3-C18 AQ

Part number: B0203F184PS030  
Dimension: 25 × 3 mm  
Eluent: ADAM-B-Phase  
Flow: 0.95 ml/min  
Detection: Coulometric  
Temperature: 28 °C  
Injection: 1 µl  
Sample:  
1: Norepinephrine  
2: Epinephrine  
3: Dihydroxybenzylamine  
4: Dopamine

